# **SPORT-ELEC®**

## **I N S T I T U T**

FR

ES

NL

GB

DE

RO

## **CEINTURE VIBRANTE PRO**

Mode d'emploi

**SPORT-ELEC® INSTITUTE**

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## WARNINGS

- ▶ Don't use while pregnant or during menstrual cycle.
- ▶ If receiving specialized treatment of any kind, please consult your doctor before attempting to use the product.
- ▶ Don't use immediately if unusual symptoms are experienced.
- ▶ Don't dismantle, modify or attempt to repair the product.
- ▶ Never use a damaged power cord, plug or power outlet.
- ▶ Avoid using directly on skin.
- ▶ Not to be used by children.
- ▶ Not to be used on the neck or joints.
- ▶ If pain develops, contact your doctor.
- ▶ Vibrating Belt is a device that utilizes effective oscillating movements in order to shape different parts of the body.
- ▶ Relieving stress and tensions. Shapes and treats all major muscle groups what leads to improved sports agility. Stimulates and massages deep tissue which aids to improve mobility.

## THE SET INCLUDES

### picture 1, page 1

- ▶ A = Vibrating Belt
- ▶ B = Remote controller
- ▶ C = AC Adaptor
- ▶ D = Connector AC Adaptor.
- ▶ 1 belt extension
- ▶ 1 carry bag

### HAND HELD CONTROL picture 2 - I, page 1

- ▶ 1 = Power button "RUN"
- ▶ 2 = "MODE" button to program each parameters
- ▶ 3 = "UP" Power increase button
- ▶ 4 = "DOWN" Power decrease button
- ▶ 5 = "HEAT" button to turn on or off the Heat function.

## FEATURES

- ▶ 4 program settings
- ▶ 5 Levels of intensity
- ▶ Heat function (between 31°C and 42°C)
- ▶ Vibrating motion
- ▶ Auto timer can be programmed from 1 to 20 mn
- ▶ Built-in auto timer will automatically shut off after program is completed, maximum 20 minutes (10 minutes by default when starting the device).

## OPERATION INSTRUCTIONS

1. Wrap the belt tightly around the body, adjust the size and fasten using the Velcro strap.
2. Attach AC adaptor to belt power connection point and plug the adaptor into the power outlet.
3. First you choose the time setting (up to 20minutes) by pressing the «UP» or «DOWN» key.
4. Then press the «MODE» button to select one of the 4 programs:
  - Program 1 = continuous vibrations
  - Program 2 = vibrations with progressive intensities
  - Program 3 = alternating vibrations (left and right) and progressive intensities

- Program 4 = intermittent vibrations with progressive intensities.

5. Once you have selected your program, press the «RUN» button to start the session. Press the «UP» or «DOWN» key to adjust intensity level, from 1 to 5.
6. Press the «HEAT» button to activate the heat. Press the button again to stop the heat. While in use the unit will heat up between 31°C and 42°C the temperature will increase once every minute. Two ventilation holes placed on top of the belt will prevent it from overheating. The unit may be stopped / powered off at any time by pressing the «RUN» key. To change the settings, the Mode button may be pressed at anytime.

## IMPORTANT PRECAUTIONS

- ▶ During normal use both the belt and adaptor will heat up.
- ▶ Unplug the unit from the power socket when not in use. Do not leave on.
- ▶ Always store in the storage bag provided.
- ▶ Never store in a damp or excessively warm environment.
- ▶ Do not put under strain or place heavy articles on top of the unit and do not sit on it.
- ▶ Clean dust or smears from the belt by wiping it with a cloth and mild detergent when it is not in use and completely cooled down.
- ▶ The Vibrating Belt will automatically switch off when it overheats. Allow it time to cool down before attempting to use again.
- ▶ Use the belt only with supplied AC adaptor
- ▶ Carefully examine the product for indications of wear and/or damage before initial use. If any such indications are found please return to supplier.

## SPECIFICATIONS

Power: Input AC100/240 V - Output DC 24V/2A  
Frequency: 50Hz - 60Hz  
Power consumption: 48W

Please apply regulations enforced in your country when disposing your product at its end of life. Therefore, please discard it in a location planned to this effect, ensuring it is disposed in an environment-friendly manner.

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## GB WARRANTY

This appliance is guaranteed for 24 months from the date of purchase. During the period of the guarantee, the distributor will take responsibility for defects occur, the appliance will be repaired or exchanged as the distributor deems necessary. Damage arising from improper use or normal wear and defects, which do not materially affect the performance of the appliance, is not covered under this guarantee. The guarantee is not valid if unauthorized repairs have been attempted or if replacement parts used are not approved by the supplier. The guarantee is valid only if it is stamped and signed by the retailer at the time of the purchase. For all claims made during the period of the guarantee, send the appliance and the warranty card to your retailer. All other claims are excluded from the guarantee unless specifically prohibited by law.

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